

# ENT Instruments Catalogue



Manufacturer & Exporters of Surgical Instruments

**darleys**  
S U R G I C A L C O.

Quality

Instruments

Since 1992



[www.darleys-surgical.com](http://www.darleys-surgical.com)

# darleys

S U R G I C A L C O.

Quality Instruments Since 1992

Manufacturer & Exporters of Surgical Instruments

2-A Kutchery Road, Sialkot-51310 Pakistan

Tel: +92 52 426 9230 Fax: +92 52 458 6215

E-Mail: [info@darleys-surgical.com](mailto:info@darleys-surgical.com)

[www.darleys-surgical.com](http://www.darleys-surgical.com)

 [darleys-surgical](https://www.linkedin.com/company/darleys-surgical)

 [fb.com/darleysurgical](https://www.facebook.com/darleysurgical)

 [twitter.com/darleysurgical](https://twitter.com/darleysurgical)

FDA cGMP CE





## Darleys Surgical Co.

Manufacturers & Exporters of  
Surgical Instruments

### Mission Statement

We, at Darleys Surgical Company strive to be a leading global manufacturer and exporters of a broad range of quality-surgical-instruments. We endeavor to provide courteous and dedicated customer service and support while supplying high-quality surgical instruments, serving a variety of disciplines, at a reasonable cost.

There is one statement that is repeated over and again by Surgeons from across the globe when referencing Darleys Instrumentation:

“Darleys” is the best choice for excellent instruments that provide safety and the best chance for procedural success.”

### COMPANY HISTORY

Since our inception in 1992, Darleys Surgical Co. has been recognized, globally, as a unique leader of innovation and design trends within the medical device industry. We have built an excellent reputation in the field by distributing a far-reaching range of surgical instruments and supplies that serve a variety of fields and disciplines.

All Darleys Instruments are masterfully hand-crafted by a team of highly-skilled and experienced craftsmen.

**All of our manufacturing facilities are ISO (9001:2008, 13485:2003) FDA cGMP and CE-Mark Certified.**

All of our factories, production facilities and distribution channels are wholly owned and operated by family members, eliminating the possibility of supply-chain issues or availability.

We invite you, our valued customer, to take the benefit from our superior instruments and services which, by themselves, reflect the premium quality and excellent value inherent to the Darleys Surgical name. It is our hope that you use this catalog as the primary purchasing source for all your Instrumentation needs. We extremely value our customer relations and greatly appreciate your business and loyalty to Darleys. We will tirelessly work to nurture and foster these relationships to create a life-long and loyal customer.

Thank You  
Darleys Surgical Co

## LIFE TIME GUARANTEE



### *Describe our Instruments*

- ✓ *Unsurpassed*
- ✓ *Ultra*
- ✓ *Reliable*
- ✓ *Great*
- ✓ *Durable*
- ✓ *Competent*
- ✓ *High-Quality*
- ✓ *Long lasting*
- ✓ *Dependable*
- ✓ *Excellent*
- ✓ *Sensational*
- ✓ *Precise*
- ✓ *Sleek*
- ✓ *Innovative*

All Instruments are being produced under the supervision of highly experienced and skilled craftsmen.

### **"WE UNCONDITIONAL GUARANTEE"**

- ❖ Use of best steel of proper specification in the production.
- ❖ Adequate testing in hardening/tempering.
- ❖ Accurate adjustment of ratchets/serrations of forceps, etc.
- ❖ Electrolytically polished, passivated and ultrasonically cleaned with best suitable chemicals.
- ❖ Flawless workmanship and superb finish.
- ❖ Minute inspection at every quality controlling process as per International Standard Organization specifications.
- ❖ Immediate settlement of claim against substandard quality.


This "GUARANTEE" is always maintained as all instruments are fabricated in uniform standard of samples exactly approved by the buyers.

Please ask for our MSDS (Material Safety Data Sheet) for better understanding on Best, Safe and long term usage of Darleys Instruments.

All Darleys Instruments are available in a high-polish or satin/matte finish according to individual preferences. Constructed from high-quality 300, 410, or 420 Series Stainless Steel, these devices meet or exceed all ASTM/AISI International Standards for Medical Device production. Undergoing the latest finishing techniques to prevent rusting and future corrosion, our state-of-the-art method includes, but is not limited to: surface grinding, broaching, heat treating, drilling, tapping, electronically passivating and I.D./O.D. Threading.

# 1 Ambulantes Ohr-Instrumentarium

## Diagnostic Ear Instruments



Gehörgangs- und Nasenspekulum  
Ear canal and nasal speculum

Art.-Nr./No. 20-100-01



Brünings Ohrtrichter  
Brünings ear speculum

Art.-Nr./No.	Größe/Size
10-183-10	0
10-183-01	1
10-183-02	2
10-183-03	3
10-183-04	4




Hartmann Ohrtrichter  
Hartmann ear speculum

Art.-Nr./No.	Größe/Size
10-200-00	<del>3,5</del> mm / 2,0 mm (Gr./
10-200-01	<del>4,5</del> mm / 3,0 mm (Gr./
10-200-02	<del>5,5</del> mm / 4,0 mm (Gr./
10-200-03	<del>6,5</del> mm / 5,0 mm (Gr./
10-200-04	<del>7,5</del> mm / 6,0 mm (Gr./



Hartmann Ohrtrichter zur Nachbehandlung, matt  
Hartmann ear speculum for post treatment, satin finished

Art.-Nr./No.	AD/OD	ID/ID
10-201-01	9 mm	8 mm
10-201-02	10 mm	9 mm



Metalstander mit 12 Aufnahmen für Ohrtrichter  
Metal stand with 12 seatings for ear speculum

Art.-Nr./No. 10-204-99

# 1 Ambulantes Ohr-Instrumentarium

## Diagnostic Ear Instruments



Farrior Ohrtrichter, rund  
Farrior ear speculum, round

Art.-Nr./No.	Innen-ø/Inner ø	Außen-ø/Outer ø
10-186-01	3,3 mm	4,3 mm
10-186-02	4,3 mm	5,1 mm
10-186-03	4,8 mm	5,6 mm
10-186-04	5,3 mm	6,1 mm
10-186-05	5,8 mm	6,6 mm
10-186-06	6,3 mm	7,1 mm
10-186-07	6,7 mm	7,5 mm
10-186-08	7,2 mm	8,0 mm
10-186-09	7,7 mm	8,5 mm



Farrior Ohrtrichter, oval, schwarz, 20°  
Farrior ear speculum, oval, black, 20°

Art.-Nr./No.	Innen-ø/Inner ø	Außen-ø/Outer ø
10-188-01	3,0 x 4,0 mm	4,0 x 5,0 mm
10-188-02	4,0 x 5,0 mm	5,0 x 6,0 mm
10-188-03	4,5 x 5,5 mm	5,5 x 6,5 mm
10-188-04	5,0 x 6,0 mm	6,0 x 7,0 mm
10-188-05	5,5 x 6,5 mm	6,5 x 7,5 mm
10-188-06	6,0 x 7,0 mm	7,0 x 8,0 mm
10-188-07	6,5 x 7,5 mm	7,5 x 8,5 mm
10-188-08	7,0 x 8,0 mm	8,0 x 9,0 mm
10-188-09	7,5 x 8,5 mm	8,5 x 9,5 mm



Farrior Ohrtrichter, oval, schwarz, 35°  
Farrior ear speculum, oval, black, 35°

Art.-Nr./No.	Innen-ø/Inner ø	Außen-ø/Outer ø
10-190-01	3,0 x 4,0 mm	4,0 x 5,0 mm
10-190-02	4,0 x 5,0 mm	5,0 x 6,0 mm
10-190-03	4,5 x 5,5 mm	5,5 x 6,5 mm
10-190-04	5,0 x 6,0 mm	6,0 x 7,0 mm
10-190-05	5,5 x 6,5 mm	6,5 x 7,5 mm
10-190-06	6,0 x 7,0 mm	7,0 x 8,0 mm
10-190-07	6,5 x 7,5 mm	7,5 x 8,5 mm
10-190-08	7,0 x 8,0 mm	8,0 x 9,0 mm
10-190-09	7,5 x 8,5 mm	8,5 x 9,5 mm




Zöllner Ohrtrichter, Leichtmetall, schwarz eloxiert  
Zöllner ear speculum, light metal, black anodized

Art.-Nr./No.	Innen-ø/Inner ø	Außen-ø/Outer ø
10-192-01	4,5 x 5,5 mm	5,5 x 6,5 mm
10-192-02	5,5 x 6,5 mm	6,5 x 7,5 mm
10-192-03	6,5 x 7,5 mm	7,5 x 8,5 mm
10-192-04	7,5 x 8,5 mm	8,5 x 9,5 mm
10-192-05	8,5 x 9,5 mm	9,5 x 10,5 mm
10-192-06	9,5 x 10,5 mm	10,5 x 11,5 mm

# 1 Ambulantes Ohr-Instrumentarium


## Diagnostic Ear Instruments



Lucae Ohrpinzette, bajonett, 14 cm  
Lucae dressing forceps, bayonet, 14 cm

Art.-Nr./No. Beschreibung/Description

- 10-272-01 **fein**, 2.0 mm
- 10-272-02 **extrafein**, 1.0 mm



Lucae Ohrpinzette, bajonett, fein,  
1:2 Zähne, 2,0 mm, 14 cm  
Lucae dressing forceps, bayonet, fine,  
1:2 teeth, 2.0 mm, 14 cm


Art.-Nr./No. 10-273-01



Troeltsch Ohrpinzette, kniegebogen, 2,0 mm  
Troeltsch dressing forceps, angled, 2.0 mm


Art.-Nr./No. Länge/length

- 10-274-00 13 cm
- 10-274-01 11 cm
- 10-274-02 15 cm
- 10-274-03 17 cm



Hartmann Ohrzängchen, gerade, gerieft, 8 x 1,8 mm, 8 cm  
Hartmann ear forceps, straight, serrated, 8 x 1.8 mm, 8 cm

Art.-Nr./No. 10-276-00




Duckbill Ohrzängchen, schnabelförmig, gerieft, 8 cm  
Duckbill ear forceps, flat jaws, serrated, 8 cm

Art.-Nr./No. 10-277-00


# 1 Ambulantes Ohr-Instrumentarium

## Diagnostic Ear Instruments




Gehörgangskürette, scharf, 17 cm  
Ear canal curette, sharp, 17 cm

Art.-Nr./No.	Größe/Size
10-669-01	1,8 mm
10-669-02	2,1 mm




Jobson Watteträger / Kürette, 18 cm  
Jobson cotton carrier / curette, 18 cm

Art.-Nr./No. 10-208-00




Snellen Ohrschlinge, gebogen, gerieft, 16 cm  
Snellen ear loop, curved, serrated, 16 cm

Art.-Nr./No.	Größe/Size
10-265-01	1
10-265-02	2
10-265-03	3




Ohr-Tamponade-Stopfer Modell Großhadern, 18 cm  
Gauze applicator, type Großhadern 18 cm

Art.-Nr./No. 10-268-18



Hartmann Tamponadezängchen, Schaft 5 cm  
Hartmann ear tampon forceps, shaft 5 cm

Art.-Nr./No.	Beschreibung/Description
10-269-01	schmal /
10-269-02	extra schmal /




Tilley Ohrtamponadezängchen, 14 cm  
Tilley dressing and sponge holding forceps, 14 cm

Art.-Nr./No. 10-271-14




## 8 Ohrzängchen

### Ear forceps




McGee Drahtschließzängchen, 0,8 x 3,5 mm, 8 cm  
 McGee wire closing forceps, 0.8 x 3.5 mm, 8 cm

Art.-Nr./No.	Beschreibung/Description
10-711-01	gerade / straight
10-711-02	rechts gebogen / curved right
10-711-03	links gebogen / curved left




Fisch-McGee Drahtschließzängchen, gerade, 0,6 x 3,5 mm, 8 cm  
 Fisch-McGee wire closing forceps, straight, 0.6 x 3.5 mm, 8 cm

Art.-Nr./No. 10-712-00



Greven Drahtschließzängchen, gerade, 0,8 x 5,0 mm, 8 cm  
 Greven wire closing forceps, straight, 0.8 x 5.0 mm, 8 cm

Art.-Nr./No. 10-714-00




Kürzzängchen für Mittelohr-Implantate, gerade,  
 1,4 x 4,0 mm, 8 cm  
 Shortening forceps for middle ear implants, straight,  
 1.4 x 4.0 mm, 8 cm

Art.-Nr./No. 10-715-00


# 1 Ambulantes Ohr-Instrumentarium

## Diagnostic Ear Instruments




Hartmann Ohrzängchen, runde Löffel, 8 cm  
Hartmann ear forceps, round cups, 8 cm

Art.-Nr./No.	Größe/Size
10-278-01	1 (1,5 mm)
10-278-02	2 (2,0 mm)
10-278-03	3 (2,5 mm)
10-278-04	4 (3,0 mm)




Hoffmann Ohrzängchen, durchschneidend, rund, 8 cm  
Hoffmann ear forceps, throughcutting, round, 8 cm

Art.-Nr./No.	Beschreibung/Description
10-279-00	2,5 mm
10-279-01	2,0 mm, gerade /
10-279-45	2,0 mm, 45° aufgebogen /



Hoffmann Ohrzängchen, durchschneidend, oval,  
2,5 x 3,5 mm, 8 cm  
Hoffmann ear forceps, throughcutting, oval,  
2.5 x 3.5 mm, 8 cm

Art.-Nr./No.	
10-280-00	




Strümpel Ohrzängchen, oval gefenstert, 5,0 x 2,5 mm, 8 cm  
Strümpel ear forceps, oval fenestrated, 5.0 x 2.5 mm, 8 cm

Art.-Nr./No.	Winkel/Angle
10-281-00	0°
10-281-45	45°
10-281-90	90°


# 1 Ambulantes Ohr-Instrumentarium

## Diagnostic Ear Instruments



Saugrohr, kniegebogen, zylindrisch, 9 cm  
 Suction tube, angled, cylindrical, 9 cm

Art.-Nr./No.	Ø ID	Länge/Length
10-300-10	1,0 mm	9 cm
10-300-15	1,5 mm	9 cm
10-300-20	2,0 mm	9 cm
10-300-25	2,5 mm	11 cm
10-300-30	3,0 mm	11 cm
10-300-40	4,0 mm	11 cm
10-300-50	5,0 mm	11 cm



Saugrohr, kniegebogen, 11,5 cm, konisch  
 Suction tube, angled, 11.5 cm, conical

Art.-Nr./No.	Ø ID	Länge/Length
10-302-30	3,0 mm	11,5 cm
10-302-35	3,5 mm	11,5 cm
10-302-45	4,5 mm	13 cm



Saugrohr, gerade, biegsam, konisch, 3,0 mm, 17,5 cm  
 Suction tube, straight, bendable, conical, 3.0 mm, 17.5 cm

Art.-Nr./No. 10-303-18



Wollstein Saughandgriff mit Unterbrecher, 11 cm  
 Woolstein suction handle with suction control, 11 cm

Art.-Nr./No. 10-305-11




House Saughandgriff mit Unterbrecher, 5 cm  
 House suction handle with suction control, 5 cm

Art.-Nr./No. 10-306-05


# 1 Ambulantes Ohr-Instrumentarium

## Diagnostic Ear Instruments



Zöllner Saugrohr, gewinkelt, 2,5 mm, 15 cm  
 Zöllner suction tube, angled, 2.5 mm, 15 cm

Art.-Nr./No. 10-310-25



Zöllner Aufsatzsauger mit Mandrin, 0,7 mm, 4 cm  
 Zöllner suction tip with stylet, 0.7 mm, 4 cm

Art.-Nr./No. 10-310-07



House Nadelsauger, Luer-Lock  
 House suction cannula, Luer-Lock

Art.-Nr./No.	Ø ID	Länge/Length
10-308-06	0,6 mm	7 cm
10-308-08	0,8 mm	7 cm
10-308-10	1,0 mm	7 cm
10-308-12	1,2 mm	7 cm
10-308-15	1,5 mm	7 cm
10-308-18	1,8 mm	7 cm
10-308-20	2,0 mm	7 cm
10-308-22	2,2 mm	7 cm
10-308-25	2,5 mm	7 cm
10-308-30	3,0 mm	7 cm
10-309-06	0,6 mm	6 cm
10-309-08	0,8 mm	6 cm
10-309-10	1,0 mm	6 cm
10-309-13	1,3 mm	6 cm



Schuknecht Saugrohr mit Unterbrecher, 45°, 7,5 cm  
 Schuknecht suction tube with suction control, 45°, 7.5 cm

Art.-Nr./No.	Ø ID
10-314-20	0,8 mm
10-314-22	0,7 mm
10-314-24	0,5 mm

Diagnostic Ear Instruments

Frazier suction tube, with suction control




É Ø Í Þ Ö Ñ Æ	Working Point	Working length
30-103-10	☞	ĈĈĈÖ Ö
30-103-11	☞	ĈĈĈÖ Ö
30-103-12	☞	ĈĈĈÖ Ö
30-103-13	☞	ĈĈĈÖ Ö
30-103-14	☞	ĈĈĈÖ Ö
30-103-15	☞	ĈĈĈÖ Ö
30-103-16	☞	ĈĈĈÖ Ö
30-103-17	☞	ĈĈĈÖ Ö
30-103-18	☞	ĈĈĈÖ Ö
30-103-19	☞	ĈĈĈÖ Ö
30-103-20	☞	ĈĈĈÖ Ö
30-103-21	☞	ĈĈĈÖ Ö
30-103-22	☞	ĈĈĈÖ Ö
30-103-23	☞	ĈĈĈÖ Ö
30-103-24	☞	ĈĈĈÖ Ö
30-103-25	☞	ĈĈĈÖ Ö
30-103-26	☞	ĈĈĈÖ Ö
30-103-27	☞	ĈĈĈÖ Ö
30-103-28	☞	ĈĈĈÖ Ö
30-103-29	☞	ĈĈĈÖ Ö
30-103-30	☞	ĈĈĈÖ Ö
30-103-31	☞	ĈĈĈÖ Ö
30-103-32	☞	ĈĈĈÖ Ö
30-103-33	☞	ĈĈĈÖ Ö
30-103-34	☞	ĈĈĈÖ Ö
30-103-35	☞	ĈĈĈÖ Ö
30-103-36	☞	ĈĈĈÖ Ö
30-103-37	☞	ĈĈĈÖ Ö
30-103-38	☞	ĈĈĈÖ Ö
30-103-39	☞	ĈĈĈÖ Ö
30-103-40	☞	ĈĈĈÖ Ö
30-103-41	☞	ĈĈĈÖ Ö
30-103-42	☞	ĈĈĈÖ Ö
30-103-43	☞	ĈĈĈÖ Ö
30-103-44	☞	ĈĈĈÖ Ö
30-103-45	☞	ĈĈĈÖ Ö
30-103-46	☞	ĈĈĈÖ Ö
30-103-47	☞	ĈĈĈÖ Ö

# 1 Ambulantes Ohr-Instrumentarium

# Otological surgery

## Diagnostic Ear Instruments



Verhoeven Saugrohr mit Griffplatte, 9 cm, gerade  
Verhoeven suction tube with finger plate, 9 cm, straight


Art.-Nr./No.	Ø ID
10-313-12	1,2 mm ○
10-314-13	1,5 mm ○
10-313-20	2,0 mm ○



Plester Saugrohr mit Griffplatte und Unterbrecher, 13 cm  
Plester suction tube with suction control, 13 cm

Art.-Nr./No.	Ø ID
10-315-15	1,5 mm ○
10-315-20	2,0 mm ○
10-315-25	2,5 mm ○
10-315-30	3,0 mm ○

## Baron suction tube with suction control



Art.-Nr./No.	working point size	working length size
10-312-10	0.4mm	60mm
10-312-11	0.4mm	70mm
10-312-12	0.6mm	70mm
10-312-13	0.7mm	70mm
10-312-14	0.8mm	70mm
10-312-15	1.0mm	70mm
10-312-16	1.4mm	70mm
10-312-17	1.6mm	70mm
10-312-18	2.0mm	70mm
10-312-19	0.4mm	80mm
10-312-20	0.7mm	80mm
10-312-21	1.4mm	80mm
10-312-22	3.0mm	80mm
10-312-23	3.0mm	80mm
10-312-24	2.0mm	100mm

# darleys

**S U R G I C A L C O .**

Quality Instruments Since 1992

Manufacturer & Exporters of Surgical Instruments

2-A Kutchery Road, Sialkot-51310 Pakistan

Tel: +92 52 426 9230 Fax: +92 52 458 6215

E-Mail: [info@darleys-surgical.com](mailto:info@darleys-surgical.com)

 [darleys-surgical](https://www.linkedin.com/company/darleys-surgical)

 [fb.com/darleyssurgical](https://www.facebook.com/darleyssurgical)

 [twitter.com/darleyssurgical](https://twitter.com/darleyssurgical)

[www.darleys-surgical.com](http://www.darleys-surgical.com)



# ENT Instruments Catalogue



Manufacturer & Exporters of Surgical Instruments

**darleys**  
SURGICAL CO.

Quality

Instruments

Since 1992

[www.darleys-surgical.com](http://www.darleys-surgical.com)