



Page
01–09

○ Micro Surgery Instruments

VESSEL DILATORS

Vessel Dilator Angled

- 11cm/4¹/₄" length
 - Tips 0.2mm Ø
- MM101-11**




Vessel Dilator Angled

- 11cm/4¹/₄" length
 - Tips 0.3mm Ø
- MM102-11**



Vessel Dilator Straight

- 13.5cm/5³/₈" length
 - Tips 0.2mm Ø
- MM103-13**



Vessel Dilator Straight

- 13.5cm/5³/₈" length
 - Curved Tips 0.2mm Ø
- MM104-18**



Vessel Dilator Long

- 18cm/7" length
 - Tips 0.3mm Ø
- MM105-18**



Vessel Dilator Long

- 18cm/7" length
 - Curved Tips 0.2mm Ø
- MM106-18**



Round Bodied Vessel Dilator

- 15cm/6" length
- Straight tips 0.3mm Ø



- Balanced Forceps



MM107-15

Round Bodied Vessel Dilator

- 15cm/6" length
- Curved tips 0.2mm Ø



- Balanced Forceps




MM108-15

Round Bodied Vessel Dilator

- 18cm/7" length
- Straight tips 0.3mm Ø



- Balanced Forceps



MM109-18

Round Bodied Vessel Dilator

- 18cm/7" length
- Curved tips 0.3mm Ø
- Balanced Forceps



MM110-18

MICRO INSTRUMENTS

Jewellers Forceps

- 11cm/4¹/₄" length
- Tips 0.3mm

MM111-11



Jewellers Forceps Long

- 15cm/6" length
- Tips 0.3mm

MM112-15



Jewellers Forceps

- 12cm/4³/₄" length
- Tips 0.3mm

MM113-12



Jewellers Forceps

- 13.5cm/5³/₈" length
- Tips 0.3mm

MM114-13



Pierse Forceps

- 12cm/4³/₄" length
- Tips 0.3mm

MM115-12



Pierse Forceps have a micro hole in the tips which securely grip tissue with minimal trauma

Micro Adson Fixation Forceps

- 12cm/4³/₄" length
- 1 x 2 Micro Teeth
- Tips 0.3mm

MM116-12



Round Bodied Forceps

- 15cm/6" length
- Straight, with platform
- Tips 0.3mm Ø
- Balanced Forceps



MM117-15

Round Bodied Forceps

- 15cm/6" length
- Curved, with platform
- Tips 0.3mm Ø
- Balanced Forceps



MM118-15

Round Bodied Angled Forceps

- 15cm/6" length
- Angled tips 0.2mm Ø
- Balanced Forceps



MM119-15

Round Bodied Forceps

- 18cm/7" length
- Straight, with platform
- Tips 0.5mm Ø
- Balanced Forceps



MM120-18

Round Bodied Forceps

- 18cm/7" length
- Curved, with platform
- Tips 0.5mm Ø
- Balanced Forceps



MM121-18

Round Bodied Angled Forceps

- 18cm/7" length
- Angled tips 0.2mm Ø
- Balanced Forceps



MM122-18


MICRO SCISSORS

Straight Scissors Flat Handle

- 15cm/6" length
- Blades 19mm long

MM123-15

Adventitia Scissors with pointed tips



MM123-15S


Available with Micro Serrations for nerve cutting - **Straight**

Curved Scissors Flat Handle

- 15cm/6" length
- Blades 19mm long

MM124-15

Dissecting Scissors with rounded tips



MM124-15S


Available with Micro Serrations for nerve cutting - **Curved**

Straight Scissors Flat Handle

- 15cm/6" length
- Blades 15mm long

MM125-15

Adventitia Scissors with pointed tips



Curved Scissors Flat Handle

- 15cm/6" length
- Blades 15mm long

MM126-15

Dissecting Scissors with rounded tips



Straight Scissors Round Handle

15cm/6" **MM127-15**


- Blades 12mm long

18cm/7" **MM128-18**

- Blades 14mm long



Adventitia Scissors with pointed tips



Curved Scissors Round Handle

15cm/6" **MM129-15**

- Blades 12mm long

18cm/7" **MM130-18**

- Blades 14mm long



Dissecting Scissors with rounded tips



Blondeel Scissors Round Handle

- Blades 15mm long
- Curved blades

15cm/6" **MM130-15**


18cm/7" **MM131-18**



MICRO NEEDLE HOLDERS

Barraquer Micro Needle Holder

- Round Bodied
- 0.4mm Jaws



Without Lock	Straight	Curved
13cm/5"	MM132-13	MM135-13
15cm/6"	MM133-15	MM136-15
18cm/7"	MM134-18	MM137-18

Castroviejo Micro Needle Holder

- Flat Handled
- 0.4mm Jaws



Without Lock	Curved
14cm/5½"	MM138-14
18cm/7"	MM139-18

Double Action Micro Needle Holder

- Round Bodied
- 0.4mm Jaws



Without Lock	Curved
15cm/6"	MM140-15
18cm/7"	MM141-18

- Sterile, single use, disposable
- Atraumatic jaws
- Colour coded for better differentiation
- Consistent and reliable closing pressure
- Designed for use on veins and arteries
- Light weight
- Blue, black and purple finish provides excellent contrast in the field
- Radiopaque
- Eliminates cleaning and processing problems
- Economical


Single Clamps
packed in boxes of 24



Double Clamps
packed in boxes of 12




MICRO CLAMPS



“A” or “V” Pattern?

Throughout the series, the “V” pattern is the all-purpose clamp, commonly used on all veins, and almost all arteries. The “A” pattern, available in all sizes, is designed for exceptionally thick-walled vessels and slippery arteries. The slight incurve of the tips of the “A” pattern jaw provides extra security against slippage.

- Light and compact design
- Extremely durable
- A unique gripping surface
- Smooth sliding bar action on all approximator clamps
- Clamp profile with parallel sided jaws
- Matt or black finish



Ordering details

SINGLE CLAMPS

MM148-396 V CLAMP SIZE
MM149-397 A B-1
0.4 - 1.0mm

MM150-398 V B-2
MM151-399 A 0.6 - 1.5mm

MM152-400 V B-3
MM153-401 A 1.0 - 2.25mm

MM154-325 HD-S
1.5 - 3.5mm

MM155-286 RD-S
2.0 - 5.0mm

DOUBLE APPROXIMATING CLAMPS

MM156-414 V CLAMP SIZE
MM157-415 A ABB - 11
0.4 - 1.0mm

MM158-416 V ABB - 2 2
MM159-417 A 0.6 - 1.5mm

MM160-418 V ABB - 3 3
MM161-419 A 1.0 - 2.25mm

MM162-326 HD - D
1.5 - 3.5mm

MM163-287 RD - D
2.0 - 5.0mm

DOUBLE CLAMPS WITH FRAME

MM164-408 V CLAMP SIZE
MM165-409 A ABB - 1
0.4 - 1.0mm

MM166-410 V ABB - 2
MM167-411 A 0.6 - 1.5mm

MM168-412 V ABB - 3
MM169-413 A 1.0 - 2.25mm

