


Page 01-07


- Plaster Instruments



TU-001
ANGLE
130 mm

TU-002
DUROTIP
ANGLE
130 mm


for soft wire up to diam.	0,5 mm	0,6 mm
for hard wire up to diam.	0,3 mm	0,4 mm



TU-003
120 mm


TU-004
DUROTIP
120 mm

for soft wire up to diam.	1,2 mm	1,2 mm
for hard wire up to diam.	0,8 mm	0,8 mm



TU-005
HARVEY
125 mm

0,4 mm
0,3 mm



TU-006
150 mm

for soft wire up to diam.	1,5 mm
for hard wire up to diam.	1,0 mm



TU-007
160 mm

1,5 mm
1,2 mm



TU-008
140 mm

1,5 mm
1,0 mm



WIRE CUTTING PLIERS, LATERAL CUTTING ACTION



TU-009
180 mm

for soft wire
up to diam. 2,0 mm
For hard wire 1,5 mm
up to diam
specially hardened
cutting edges, front
and lateral cutting
action



TU-010
REILL
175 mm

for soft wire
up to diam. 2,0 mm
For hard wire 1,6 mm
up to diam
specially hardened
cutting edges, for
cutting off wires close
to the bone



TU-011
235 mm

for soft wire
up to diam. 3,0 mm
For hard wire 2,2 mm
up to diam



KILLIAN
TU-012
200 mm



TU-013 - 014
TU-013 195 mm
TU-014 225 mm



LEBSCHE
TU-015
Sternum chisel
200 mm



TU-016
230 mm




COLLIN
TU-017
195 mm



GLUCK
TU-018
190 mm



PLASTER SHEARS



TU-019
STILLE-AESULAP
260 mm



TU-020
STILLE-AESULAP
360 mm

Blade with reinforced back



TU-021
STILLE
230 mm

TU-022
Spare blade



TU-023
STILLE
370 mm

TU-024
Spare blade



TU-025
STILLE
240 mm

TU-026
Spare blade



PLASTER SPREADER, PLASTER BREAKERS



BEESON
TU-027
30 cm



MOD.USA
TU-028
23 cm



HAGLUND-STILLE
TU-029
20 cm



TU-030
29 cm



WOLFF
TU-031
10 cm



WOLFF
TU-032
24 cm