

Page 01-42

- Plastic Surgery Instruments

Gorney Face Lift Scissors

Tungsten Carbide Edges

Straight	Curved	Length
MPL101-17	MPL104-17	17cm/6 ³ / ₄ "
MPL102-19	MPL105-19	19cm/7 ¹ / ₂ "
MPL103-23	MPL106-23	23cm/9"

Flat tips



Gorney-Freeman Face Lift Scissors

Tungsten Carbide Edges

Open shanks

Straight	Curved	Length
MPL107-17	MPL110-17	17cm/6 ³ / ₄ "
MPL108-19	MPL111-19	19cm/7 ¹ / ₂ "
MPL109-23	MPL112-23	23cm/9"

Flat tips



Kaye Face Lift Scissors

Tungsten Carbide Edges

Curved

Code	Length
MPL113-17	17cm/6 ³ / ₄ "
MPL114-19	19cm/7 ¹ / ₂ "
MPL115-23	23cm/9"



Aufricht Scissors

Tungsten Carbide Edges

Curved

14cm/5¹/₄"

MPL116-14



Goldman-Fox Scissors

Tungsten Carbide Edges

Tapered blades to fine tips

Curved

13cm/5"

MPL117-13



DISSECTING SCISSORS

Ragnell (Kilner) Scissors

Tungsten Carbide Edges
Curved flat tips

Code	Length
MPL118-13	13cm/5 ¹ / ₄ "
MPL119-15	15cm/6"
MPL120-18	18cm/7"



Greenberg PAR Scissors

Tungsten Carbide Edges
Curved tapered flat tips
11.5cm/4³/₄"

MPL121-11



Littler Plastic Surgery Scissors

Curved
12cm/4³/₄"
Hole in the blades for the suture to be threaded through

MPL122-12



Kahn Face Lift Scissors

Super Cut Edges
Curved
14.5cm/5³/₄"

MPL123-14



Fomon Saber Face Lift Scissors

Curved, round blunt tips
Super Cut Edges
13.5cm/5¹/₄"

MPL124-13



Stevens Tenotomy Scissors

11.5cm/4½" Curved

MPL125-11

Jameson-Werber Scissors

13cm/5¼" Curved

MPL126-13

Jameson Scissors

15cm/6" Curved

MPL127-15

Jameson Scissors

15cm/6" Curved

Strong pattern

MPL128-15

Crayford-Reynolds Scissors

18cm/7" Curved

MPL129-18

METZENBAUM SCISSORS

Waldram Scissors (with larger bows)
12.5cm/5" Curved

MPL130-12



Slim Baby Metzenbaum Scissors
14cm/5½" Curved

MPL131-14



Slim Metzenbaum Scissors
18cm/7" Curved

MPL132-18



Baby Metzenbaum Scissors
14cm/5½" Curved

MPL133-14



McIndoe / Metzenbaum Scissors
18cm/7" Curved

MPL134-18



Derf Needle Holder

12cm/4³/₄" TC



MPL135-12

Smooth



MPL136-12



Halsey Needle Holder

13cm/5" TC



MPL137-13

Smooth



MPL138-13



Crile Wood Needle Holder TC



MPL139-1

5 15cm/6"



MPL140-1

8 18cm/7"



MPL141-20


20cm/8"



NEEDLE HOLDERS

Nieverts Needle Holder (offset ring)

13cm/5¹/₄" TC



MPL142-13



Foster Gillies Combined Needle Holder/Scissors


13cm/5" TC



MPL143-13


Gillies Combined Needle Holder/Scissors

16cm/6¹/₄" TC



MPL144-16


Gillies Dissecting Forceps

Toothed
15cm/6"**MPL145-15**


McIndoes Dissecting Forceps

Non Toothed
15cm/6"**MPL146-15**

Debakey-Diethrich Tissue Forceps

1.5mm Tips
15cm/6"**MPL147-15**

Debakey Tissue Forceps

2mm Tips
15cm/6"**MPL148-15**

Bipolar Forceps

1mm Tips
11cm/4¹/₂"**MPL149-11**

Bipolar Forceps

1mm Tips
16cm/6¹/₄" Length**MPL150-16**

DISSECTING FORCEPS

Adson Forceps

Non Toothed
12cm/4³/₄"

MPL151-12



Adson Forceps

Toothed
12cm/4³/₄"

MPL152-12



Micro Adsons Forceps

Non Toothed
12cm/4³/₄"

MPL153-12



Micro Adson Forceps

Toothed
12cm/4³/₄"

MPL154-12



Adson Brown Forceps

7x7 Teeth
12cm/4³/₄"

MPL155-12



Micro Fixation Forceps

Micro Teeth
12cm/4³/₄"

MPL156-12



Kilner Detachable Skin Hook

15cm/6"

MPL157-15

Gillies Skin Hook

18cm/7"

MPL158-18 Small size**MPL159-07** Large size

Klienert-Kutz Skin Hook

16cm/6¹/₄"**MPL160-01**

Joseph Double Skin Hook

16cm/6¹/₄"**MPL161-02**

Catspaw Retractors

16cm/6¹/₄"

Semi-sharp Prongs

MPL162-03

Senn-Miller Double Ended Retractor

16cm/6¹/₄"**MPL163-01** Sharp Prongs**MPL164-02** Blunt Prongs

Kilner-Miniature Pattern Catspaw Retractor

15cm/6"

MPL165-01

RETRACTORS

Ragnell-Davies Retractor

14cm/5½"

MPL166-14



Mini-Langenbeck Retractor

16cm/6¼"

MPL167-10 10x6mm

MPL168-17 17x5mm



Desmarres Saddle Retractor

16cm/6¼"

12mm width

MPL169-12



Alms Retractor

4 x 4 Teeth Blunt

7cm/2¾"

4mm depth Prongs

MPL170-70



Baby Weitlaner S/R Retractor

2 x 3 Blunt Prongs

11.5cm/4½"

MPL171-10



West's S/R Retractor

3 x 4 Sharp Prongs

14cm/5½"

MPL172-14



Barrons Scalpel Blade Handle No.3

13cm/5¹/₄"

MPL173-13



Beaver Blade Holder

13cm/5¹/₄"

MPL174-13



Eckoff Mapping Pen

14.5cm/5³/₄"

MPL175-14



Sommerlad Marking Pen

15cm/6"

MPL176-15



Stainless Steel Measuring Ruler

Code	Length
MPL177-15	15cm/6"
MPL178-30	30cm/12"



Mitchels Trimmer

16cm/6³/₄"

MPL179-16



Castroviejo Measuring Caliper

7cm/2³/₄"

Measuring range 0 to 25mm

MPL180-10



SKIN GRAFT KNIVES

Watson Skin Graft Knife

30.5cm/12"

Right handed

MPL181-30

Left handed

MPL182-30



The roller bar is smooth and does not oscillate or rotate in use. The graft thickness is controlled by slackening the locking nut and rotating the adjuster nut until the required thickness and then set by tightening the locking nut. The graft thickness is consistent as the roller bar can only be adjusted in parallel.

Braithwaite Skin Graft Knife


30.5cm/12"

Right handed

MPL183-30

Left handed

MPL184-30



The grooved roller bar oscillates from side to side as well as rotating. The graft thickness is controlled by rotating the adjuster nuts at the end of knife, this will move the roller bar up or down.

Cobbett Skin Graft Knife


30.5cm/12"

Right handed

MPL185-30

Left handed

MPL186-30



The grooved roller bar only oscillates from side to side. The adjustment for the thickness of the graft is the same as the Braithwaite knife.

Humby Skin Graft Knife


36.5cm/14¹/₄"

Right handed

MPL187-36

Left handed

MPL188-36



The grooved roller bar rotates and oscillates from side to side. The graft thickness is controlled by adjusting the thumb screws on top of the support pillars.

Silver Skin Graft Knife

19cm/7¹/₂"

Right handed

MPL189-19

Left handed

MPL190-19



Knife	Blades code
Cobbett	MPL191
Watson	MPL192
Braithwaite	MPL193
Silver	MPL194



Tufnol Skin Graft Board

200mm x 100mm / 8 x 4"
Bevelled edge

MPL195-20

Assumpcao Face Lift Forceps

For precise marking of dissected facial flaps for excision
16cm/6¹/₄"


Without ratchet

MPL196-16

With ratchet

* **MPL197-16**

illustrated with ratchet*



Pitanguy Forceps

Flap demarcator for accurate marking of facial flaps for excision

Code

Length

* **MPL198-15**

15cm/6"

MPL199-26

26cm/10¹/₄" (for abdominal flaps)

Illustrated *



Marten Flap Forceps

For accurate marking of skin flaps
14.5cm/5³/₄"

Without ratchet

* **MPL200-14**

With ratchet

MPL201-14

Illustrated *



Martel Flap Grasping Forceps

14cm/5¹/₄"

Two pins in opposing holes on the Jaws

MPL202-14



BREAST MARKERS

Freeman Areola Marker

Code	Diameter
MPL203-36	36mm Ø
MPL204-38	38mm Ø
MPL205-40	40mm Ø
MPL206-42	42mm Ø
MPL207-45	45mm Ø



Freeman Double Ended Areola Marker

Code	Diameter
MPL208-38	38/42mm Ø
MPL209-42	42/45mm Ø



Galvao Nipple Marker

Code	Diameter
MPL210-40	40mm Ø standard
MPL211-50	50mm Ø large



McKissock Keyhole Breast Reduction Marker


Code	Diameter
MPL212-38	38mm Ø
MPL213-42	42mm Ø
MPL214-45	45mm Ø



- Precisely marks amount of Tissue for removal
- Used together with Freeman Marker
- Arms are graduated in centimeters for marking tissue to be removed

Tebbetts Fibre Optic Breast Retractor

With fibre optic illumination and suction



- Fibre optic illumination with suction for smoke evacuation
- Angled Retractor tip with blunt teeth
- Polished mirror finish on blade to reflect light inside the pocket

Code	With Teeth	
	Length	Width
MPL215-260	90mm	24mm
MPL216-260	90mm	30mm
MPL217-260	150mm	36mm
MPL218-260	190mm	36mm

Code	Blunt	
	Length	Width
MPL219-260	90mm	24mm
MPL220-260	90mm	30mm
MPL221-260	150mm	36mm
MPL222-260	190mm	36mm



FIBRE OPTIC LIGHT LEAD

Fibre Optic Light Lead

Code	Length	Diameter
MPL223-180	1800mm	3.5mm Ø
MPL224-230	2300mm	3.5mm Ø



Giunta Nasal Scissors

Tungsten Carbide Edges
Double angled blades
13.5cm/5¹/₄"

MPL225-13



Daniels Nasal Scissors

40° angled blades
10cm/4"

MPL226-10

For dissection of Nasal Tip Cartilage



Foman Upper Lateral Scissors

Curved
13.5cm/5¹/₄"

MPL227-13



Foman Dorsal Scissors

Tungsten Carbide Edges
13cm/5"

MPL228-13



Cottle Dorsal Scissors

Tungsten Carbide Edges
16cm/6¹/₄"

MPL229-16





NASAL INSTRUMENTS

Knight Nasal Scissors

Tungsten Carbide Edges
18cm/7"

MPL230-18



Heymann Nasal Scissors

18cm/7"

MPL231-18



Caplan Nasal Shears

Double action angled
20cm/8"

MPL232-20



Cottle Septum Knife

14cm/5½"

MPL233-14

Joseph button end Knife

Straight

14cm/5½"

MPL234-14

Convex button end Knife

Curved

16cm/6¾"

MPL235-16

Freer Mucosa Knife

15cm/6"

MPL236-15

Joseph Angled button end Knife

17.5cm/7"

MPL237-17

Freer Elevator

19cm/7½"

MPL238-18

NASAL INSTRUMENTS

MacDonalds Elevator

19cm/7½"

MPL239-19



Watson Cheyne Elevator

19cm/7½"

MPL240-18



Cottle Septum Elevator

Double Ended graduated
20cm/8"

MPL241-20



Masing Septum Elevator

20cm/8"

MPL242-20



Carter Septum Elevator

Double Ended with curette
21cm/8¼"

MPL243-21



Killian Septum Elevator

Double Ended
19cm/7½"

MPL244-19



Kilner Alae Retractor

85mm

Code	Blunt Prongs	Code	Sharp Prongs
MPL245-10	10mm	MPL248-10	10mm
MPL246-13	13mm	MPL249-13	13mm
MPL247-15	15mm	MPL250-15	15mm



Cottle-Neivert Alae Retractor

2 prong blunt ball tips
with knife guard
20cm/8"

MPL251-20



Foman Alae Retractor

2 prong blunt ball tips
16cm/6 1/4"

MPL252-16



Cottle Hook

2 prong blunt ball and sharp tips
14cm/5 1/2"

Sharp prongs

MPL253-14



Ball on right

MPL254-14



Ball on left

MPL255-14



Cottle Retractor

15cm/6"

MPL256-15



Aufricht Nasal Retractor

16cm/6 1/4"
45mm Blade

MPL257-16



Aufricht Nasal Retractor


16cm/6 1/4"
with Fibre Optic illumination
45mm Blade

MPL258-16




NASAL INSTRUMENTS

Killian Speculum
13cm/5"
Plain blade



Code	Blade
MPL259-35	35mm
MPL260-55	55mm
MPL261-65	65mm
MPL262-75	75mm
MPL263-90	90mm

Cottle Speculum
14cm/5½"
Slender blade




Code	Blade
MPL264-20	20mm
MPL265-30	30mm
MPL266-40	40mm
MPL267-60	60mm

Joseph Saw

19cm/7 1/2"




Illustrated *



Straight
MPL268-19



Cranked to right
MPL269-19



Cranked to left
MPL270-19 *

Joseph Saw

Horizontal shaft
125mm x 85mm

Cranked right
MPL271-19



Illustrated *

Cranked left
* **MPL272-19**

Ballenger Swivel Knife

Bayonet shaped
19.5cm/7 1/2"



Code	Width
MPL273-19	3mm
MPL274-19	4mm
MPL275-19	5mm

Jansen Middleton Septum Forceps

19cm/7 1/2"
Double action


MPL276-19






NASAL INSTRUMENTS


Walsham Septum Straightening Forceps
23cm/9"



Universal
MPL277-23



Left
MPL278-23



Right
MPL279-23

Asch Septum Straightening Forceps
23cm/9"



MPL280-23

Toffee Mallet

19cm/7½"
Head size 22mm Ø
Weight 165g/6oz

MPL281-21



Ultra Fine Mallet

17cm/6¾"
Head size 20mm Ø
Weight 145g/5oz
Tufnol Handle

MPL282-17



Read Nasal Gouge

18cm/7"

Code	Width
MPL283-04	4mm
MPL284-05	5mm
MPL285-06	6mm



Killian Septum Gouge

Round cutting edge bayonet shaped 6mm width
16cm/6¼"

MPL286-06



Killian Septum Gouge

Notched edge bayonet shaped 6mm width
16cm/6¼"

MPL287-06



TNASAL INSTRUMENTS

Matthew Osteotome

18cm/7"
11mm wide
Bevelled at one corner

MPL288-11



Gubisch Osteotome

18.5cm/7¼"
3mm wide
Aluminium handle

MPL289-03



Rubin Osteotome

18cm/7"

Code	Size
MPL290-08	8mm
MPL291-10	10mm
MPL292-12	12mm
MPL293-14	14mm
MPL294-16	16mm



French (Gillies) Osteotome

14cm/5½"

Code	Size
MPL295-02	2mm
MPL296-05	5mm
MPL297-07	7mm
MPL298-09	9mm
MPL299-11	11mm
MPL300-13	13mm



Cinelli Double Guarded Osteotome

16cm/6¼"

Code	Size
MPL301-08	8mm
MPL302-10	10mm
MPL303-12	12mm
MPL304-14	14mm
MPL305-16	16mm



Jays Osteotome

18cm/7"

Code	Size
MPL306-02	2mm
MPL307-04	4mm
MPL308-05	5mm
MPL309-06	6mm
MPL310-08	8mm
MPL311-10	10mm
MPL312-12	12mm



Cottle Chisel

18.5cm/7 1/2"

Code	Size
MPL313-04	4mm
MPL314-06	6mm
MPL315-07	7mm
MPL316-08	8mm
MPL317-09	9mm
MPL318-12	12mm



McIndoe Chisel

15cm/6"

Code	Size
MPL319-09	9mm
MPL320-11	11mm
MPL321-13	13mm
MPL322-15	15mm



Sheehan Osteotomes

15cm/6"

Set of 5 sizes
2mm to 10mm

Code	Size
MPL323- 02	2mm
MPL324- 04	4mm
MPL325- 06	6mm
MPL326- 08	8mm
MPL327-10	10mm



NASAL INSTRUMENTS

Braithwaite Chisel

14.5cm/5³/₄"
5mm wide

Code

- MPL328-14** Without guard
- MPL329-14** Guard on left
- MPL330-14** Guard on right



Foman Chisel with guard

16.5cm/6¹/₂"

Code

Size

- MPL331-04** 4mm
- MPL332-05** 5mm
- MPL333-06** 6mm
- MPL334-07** 7mm



Silver Chisel with guard (Frenchay)

18cm/7"

Code

- MPL335-18** Straight
- MPL336-18** Curved left
- MPL337-18** Curved right



Nievert Chisel with Guard

20cm/8"



Code

- MPL338-20** Straight
- MPL339-20** Curved left
- MPL340-20** Curved right





Foman Double Ended Nasal Rasp

21cm/8¹/₄" Tungsten Carbide Surface





MPL341-8012

Course



MPL342-8056

Medium



MPL343-8091

Fine

Parkes Nasal Rasp

21cm/8¹/₄" Tungsten Carbide Surface
Backward cutting



MPL344-8001

Course



MPL-345-8005

Medium



MPL-346-8010

Fine

Maltz Nasal Rasp

18.5cm/7¹/₂"

MPL347-18



Maltz Nasal Rasp

18.5cm/7¹/₂"

Tungsten Carbide Surface

MPL348-18





NASAL INSTRUMENTS

Peet Nasal Rasp - Curved

17.5cm/6³/₄"

Diamond Surface

MPL349-17



Aufricht (Peet) Nasal Rasp

22cm/8³/₄"

Diamond Surface

MPL-350-22



Glabella Nasal Rasp

21cm/8¹/₄"

Tungsten Carbide Surface

MPL351-21



Foman Double Ended Nasal Rasp

21cm/8¹/₄"

Diamond Surface


MPL352-21





Millar Nasal Rasp

18cm/7"


Tungsten Carbide Surface



MPL353-18



MPL354-18



Frontotemporal Dissector22.5cm/8³/₄" Straight 10mm wide**MPL355-22**

For initial dissection through forehead incisions

Frontotemporal Dissector24cm/9¹/₂" Curved 7mm wide**MPL356-24**

For initial dissection

Periosteal Dissector25cm/9³/₄" Curved 9mm wide**MPL357-25**

For a posterior scalp dissection

Periosteal Dissector/Elevator24cm/9¹/₂" S - Shaped 9mm wide**MPL358-24**

For elevation of the periosteum

Frontoglabellar Dissector26cm/10¹/₂" S - Shaped 7mm wide**MPL359-26**

Used in the dissection of the naso-glabellar area of the forehead

BROWLIFT INSTRUMENTS

Temporal Line Dissector

23.5cm/9¹/₄" Straight 10mm wide

MPL360-23



Separates the temporal fascia

Zygomatic Arch Elevator

22.5cm/8³/₄" Straight 5mm wide

MPL361-22



For use in the lower eyelid approach

Scalp Elevator

24cm/9¹/₂" Straight 7mm wide

MPL362-24



Scalp Elevator

25.5cm/10" Curved 12mm wide

MPL363-25



Scalp Elevator

24cm/9¹/₂" Half Curved with angled tip 12mm wide

MPL364-24



Flap Dissector

22.5cm/8³/₄" Curved 6mm wide

MPL365-22



For dissection of the deep temporal fascia

Arcus Marginalis Dissector

24cm/9¹/₂" S - Shaped 8mm wide

MPL366-24



For dissection along the orbital rim

Endoforehead Elevator

24cm/9¹/₂" Half Curved 7mm wide

MPL367-24



Transoral Dissector

23.5cm/9¹/₄" 7mm wide

Curved tip

MPL368-23



Orbital Rim Dissector

24cm/9¹/₂"
Half Curved

MPL369-24





BROWLIFT INSTRUMENTS

Mini Forehead Elevator
18cm/7" Shaft
Slightly Curved

MPL370-18



Mini Forehead Elevator
18cm/7" Shaft
Half Curved

MPL371-18



Mini Forehead Dissector
18cm/7"
Full Curved

MPL372-18




Shaper Hook Sissors

2.7mm diameter
Insulated shaft

MPL373-800 Curved left

MPL374-700 Curved right




Shaper Grasping Forceps

2.7mm diameter
Insulated shaft

MPL375-800 Curved left

MPL376-900 Curved right




Shaper Punch Forceps

2.7mm diameter
Insulated shaft

MPL377-800 Curved left

MPL378-900 Curved right



BROWLIFT INSTRUMENTS

Autoclavable Telescope
4mm diameter, 30° angle

MPL379-030



Container for Telescope
600-290-052

Sheath for Telescope

For 4mm Telescope
With 2 rotating stopcocks and
flat spatula tip

MPL380-175



Sheath for Telescope

For 4mm Telescope
With one rotating stopcock and
curved spatula tip

MPL381-180



Sheath for Telescope

For 4mm Telescope
With one rotating stopcock and
angled spatula tip

MPL382-180



Vessel Dilator

11cm Angled
Tip diameter 0.3mm

MPL383-03**Vessel Dilator**

11cm Angled
Tip diameter 0.2mm

MPL384-02**Round Bodied Vessel Dilator**

15cm Straight
Tip diameter 0.3mm

MPL385-03**Round Bodied Vessel Dilator**

18cm Straight
Tip diameter 0.3mm

MPL386-03**Round Bodied Vessel Dilator**

18cm Curved
Tip diameter 0.3mm

MPL387-03

MICRO INSTRUMENTS

Jewellers Forceps

11cm
Tip diameter 0.3mm



MPL388-03

Jewellers Forceps

12cm
Tip diameter 0.3mm



MPL389-03

Pierse Forceps

12cm
Tip diameter 0.3mm



MPL390-03

Pierse Forceps have a micro hole in the tips which securely grip tissue with minimal trauma

Round Bodied Forceps

15cm Straight, with platform
Tip diameter 0.3mm



MPL391-03

Round Bodied Forceps

15cm Curved, with platform
Tip diameter 0.3mm



MPL392-03

Round Bodied Forceps

18cm Straight, with platform
Tip diameter 0.5mm



MPL393-05

Round Bodied Forceps

18cm Curved, with platform
Tip diameter 0.5mm




MPL394-05


Straight Scissors Flat Handle

15cm length
(blades 19mm long)

MPL395-15

Adventitia Scissors with pointed tips




	Available with Micro Serrations for nerve cutting - Straight
<p>SAS-15T</p>	


Curved Scissors Flat Handle

15cm length
(blades 19mm long)

MPL396-15

Dissecting Scissors with rounded tips




	Available with Micro Serrations for nerve cutting - Curved
<p>SDC-15T</p>	

Straight Scissors Flat Handle

15cm length
(blades 15mm long)

MPL397-15

Adventitia Scissors with pointed tips



Curved Scissors Flat Handle

15cm length
(blades 15mm long)

MPL398-15

Dissecting Scissors with rounded tips



MICRO INSTRUMENTS


Straight Scissors Round Handle

15cm length
(blades 12mm long)

MPL399-15



Adventitia Scissors with pointed tips




Curved Scissors Round Handle

15cm length
(blades 12mm long)

MPL400-15



Dissecting Scissors with rounded tips




Straight Scissors Round Handle

18cm length
(blades 14mm long)

MPL401-18



Adventitia Scissors with pointed tips



Curved Scissors Round Handle

18cm length
(blades 14mm long)

MPL402-18



Dissecting Scissors with rounded tips



Straight Needle Holder

15cm length without lock

MPL403-15

Curved Needle Holder

15cm length without lock

MPL404-15

Straight Needle Holder

18cm length without lock


MPL405-18

Curved Needle Holder

18cm length without lock

MPL406-18

MICRO CLAMPS



“A” or “V” Pattern?

Throughout the series, the “V” pattern is the all-purpose clamp, commonly used on all veins, and almost all arteries. The “A” pattern, available in all sizes, is designed for exceptionally thick-walled vessels and slippery arteries. The slight incurve of the tips of the “A” pattern jaw provides extra security against slippage.

- Light and compact design
- Extremely durable
- A unique gripping surface
- Smooth sliding bar action on all approximator clamps
- Clamp profile with parallel sided jaws
- Matt or black finish



Ordering details

SINGLE CLAMPS

MPL407-396 V
MPL408-397 A



CLAMP SIZE
B-1
0.4 - 1.0mm

MPL409-398 V
MPL410-399 A




B-2
0.6 - 1.5mm

MPL411-400 V
MPL412-401 A




B-3
1.0 - 2.25mm

MPL-413-325 V



HD-S
1.5 - 3.5mm


MPL414-286



RD-S
2.0 - 5.0mm

DOUBLE APPROXIMATING CLAMPS

MPL415-414 V
MPL416-415 A



CLAMP SIZE
ABB - 11
0.4 - 1.0mm

MPL417-416 V
MPL418-417 A




ABB - 22
0.6 - 1.5mm

MPL419-418 V
MPL420-419 A






ABB - 33
1.0 - 2.25mm

MPL421-326 V



HD - D
1.5 - 3.5mm

MPL422-287



RD - D
2.0 - 5.0mm

DOUBLE CLAMPS WITH FRAME

MPL423-408 V
MPL424-409 A



CLAMP SIZE
ABB-1
0.4 - 1.0mm

MPL425-410 V
MPL426-411 A



ABB - 2
0.6 - 1.5mm

MPL427-412 V
MPL428-413 A



ABB-3
1.0 - 2.25mm

